

## Product Datasheet

### (Pro3) GIP, human (orb1981131)

<b>Catalog Number</b>	orb1981131
<b>Category</b>	Proteins
<b>Description</b>	(Pro3)GIP, human is a stable, high-affinity ( $K_{i/Kd} = 0.90$ nM) full agonist of the human GIP receptor. It is a valuable research tool for investigating the role of the GIPR in metabolic diseases, enabling in vitro and in vivo studies relevant to obesity and diabetes.
<b>Target</b>	IGF-1R
<b>MW</b>	4951.53
<b>CAS Number</b>	299898-52-5
<b>Formula</b>	$C_{226}H_{338}N_{60}O_{64}S$
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)