

Product Datasheet

Exendin-4 (3-39) (orb1981156)

Catalog Number	orb1981156
Category	Proteins
Description	Exendin-4 (3-39) is a truncated, potent GLP-1 receptor agonist used in diabetes and HPA axis disorder research. It is widely applied in both in vitro receptor binding assays and in vivo studies investigating glucose metabolism and neuroendocrine functions.
Target	Glucagon Receptor
MW	3992.44
CAS Number	196109-31-6
Formula	$C_{176}H_{272}N_{46}O_{58}S$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)