

Product Datasheet

d[Cha4]-AVP TFA (orb1981523)

Catalog Number	orb1981523
Description	d[Cha4]-AVP TFA is a potent and selective agonist for the vasopressin (AVP) V1b receptor, with a high affinity (K_i of 1.2 nM). It shows greater selectivity for the V1b receptor over the V1a receptor, V2 receptor, and oxytocin receptors, demonstrating specificity in targeting the human V1b receptor [1] [2].
Storage	-20°C
Note	For research use only
Formula	$C_{52}H_{72}F_3N_{13}O_{13}S_2$
Purity	98.00%
MW	1208.33
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)