

Product Datasheet

UFP-803 TFA (orb1981625)

Catalog Number	orb1981625
Category	Small Molecules
Description	UFP-803 TFA, a potent ligand for the urotensin-II receptor (UT), exhibits minimal residual agonist activity, rendering it a valuable asset for exploring the UT system's function in physiology and pathology.
Target	GPCR
MW	1175.25
Formula	$C_{52}H_{65}F_3N_{10}O_{14}S_2$
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)