

Product Datasheet

Apelin-17(human, bovine) TFA (orb1981666)

Catalog Number	orb1981666
Category	Small Molecules
Description	Apelin-17(human, bovine) TFA is an endogenous agonist of the orphan G protein-coupled receptor APJ, demonstrating affinity for human APJ receptors on HEK 293 cells with a pIC50 value of 9.02. It also blocks the entry of certain HIV-1 and HIV-2 strains into NP2/CD4 cells that express the APJ receptor.
Target	HIV Protease,Apelin receptor
MW	2252.57
Formula	$C_{98}H_{157}F_3N_{34}O_{22}S$
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)