

Product Datasheet

Saralasin TFA (orb1981717)

Catalog Number	orb1981717
Category	Proteins
Description	Saralasin TFA ([Sar1,Ala8] Angiotensin II) is a synthetic octapeptide analog of angiotensin II, functioning as a competitive antagonist at angiotensin II receptors with a K_i value of 0.32 nM for 74% of its binding sites, while also exhibiting partial agonist activity. It is utilized in the study of renovascular and renin-dependent (angiotensinogenic) hypertension.
Target	RAAS
MW	1026.07
Formula	$C_{44}H_{66}F_3N_{13}O_{12}$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)