

Product Datasheet

TIP 39, Tuberoinfundibular Neuropeptide (orb1983854)

Catalog Number	orb1983854
Category	Proteins
Description	This synthetic peptide is a potent agonist for the human and rat parathyroid hormone 2 (PTH2) receptor, mimicking the endogenous tuberoinfundibular neuropeptide TIP39. It is a valuable research tool for studying neuroendocrine signaling, pain modulation, and cardiovascular function in both in vitro and in vivo experimental systems.
Target	AChR
MW	4504.2
CAS Number	277302-47-3
Formula	$C_{202}H_{325}N_{61}O_{54}S$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)