

Product Datasheet

Neuropeptide S (human) (TFA) (orb1983936)

Catalog Number	orb1983936
Category	Small Molecules
Description	Neuropeptide S human (TFA) is a potent endogenous neuropeptide S receptor agonist (EC50=9.4 nM).
MW	2301.52
Formula	$C_{95}H_{156}N_{31}F_3O_{30}S$
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)