

Product Datasheet

SPR741 TFA (1179330-52-9 free base) (orb1984073)

Catalog Number	orb1984073
Category	Proteins
Description	SPR741 TFA (NAB741 TFA) is a polymyxin B-derived cationic peptide that functions as a potentiator. It is used in research to enhance the activity of other antibiotics against Gram-negative bacteria in both in vitro and in vivo experimental models.
MW	1106.15
Formula	$C_{46}H_{74}F_3N_{13}O_{15}$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)