

Product Datasheet

Bexarotene D4 (orb1992936)

Catalog Number	orb1992936
Category	Small Molecules
Description	Bexarotene D4 is a deuterium-labeled Bexarotene (LGD1069). Bexarotene is a selective RXR agonist used in the treatment of cutaneous T-cell lymphoma.
Target	Retinoid Receptor, Autophagy
MW	352.5
CAS Number	2182068-00-2
Formula	C ₂₄ H ₂₈ O ₂
SMILES	[2H]c1c([2H])c(c([2H])c([2H])c1C(O)=O)C(=C)c1cc2c(cc1C)C(C)(C)CCC2(C)C
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558