

Product Datasheet

Donkey anti-Sheep IgG (H&L) - Affinity Pure, HRP Conjugate, Cross-Adsorbed (orb1993910)

Catalog Number	orb1993910
Category	Antibodies
Description	Donkey anti-Sheep IgG (H&L) - Affinity Pure, HRP Conjugate, Cross-Adsorbed
Clonality	Monoclonal
Species/Host	Donkey
Conjugation	HRP
Reactivity	Sheep
Tested applications	ELISA, IHC, WB
Cross Reactivity	Based on IEP, no reactivity is observed to: Â· non-immunoglobulin sheep serum proteins Â· IgG from human, mouse or rabbit
Antibody Type	Secondary Antibody
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)