

Product Datasheet

ARHGAP17 Peptide - middle region (orb2000811)

Description	ARHGAP17 Peptide - middle region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	This is a synthetic peptide designed for use in combination with ARHGAP17 Rabbit Polyclonal Antibody (orb589084). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
Protein Sequence	Synthetic peptide located within the following region: AGPSPHTLRRRAVKKPAPAPPKPGNPPPGHPGGQSSSGTSQHPPSLSPKPP
MW	45 kDa
Uniprot ID	Q68EM7
NCBI	NP_001006635.1
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)