

Product Datasheet

DRAXIN Peptide - C-terminal region (orb2005486)

Catalog Number	orb2005486
Category	Proteins
Description	DRAXIN Peptide - C-terminal region
Predicted Reactivity	Human
Form/Appearance	Lyophilized powder
Buffer/Preservatives	Lyophilized powder
Protein Sequence	Synthetic peptide located within the following region: TLDMALFDWTDYEDLKPDPGWPSAKKKEKHRGKLSSDGNETSPAEGEPCDH
UniProt ID	Q8NBI3
Tested applications	WB
Application notes	This is a synthetic peptide designed for use in combination with DRAXIN Rabbit Polyclonal Antibody (orb587294). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
NCBI	NP_940947
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)