

## Product Datasheet

### CAPN8 Peptide - middle region (orb2007118)

<b>Catalog Number</b>	orb2007118
<b>Category</b>	Proteins
<b>Description</b>	CAPN8 Peptide - middle region
<b>Predicted Reactivity</b>	Human
<b>Form/Appearance</b>	Lyophilized powder
<b>Buffer/Preservatives</b>	Lyophilized powder
<b>Protein Sequence</b>	KLIRLRNPWGEVEWSGAWSDDAPEWNHIDPRRKEELDKKVEDGEFWMSLS
<b>UniProt ID</b>	<b>A6NHC0</b>
<b>Tested applications</b>	WB
<b>Application notes</b>	This is a synthetic peptide designed for use in combination with CAPN8 Rabbit Polyclonal Antibody (orb586432). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
<b>Storage</b>	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
<b>Note</b>	For research use only
<b>NCBI</b>	<b>NP_001137434</b>
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)