

Product Datasheet

RAB37 Peptide - middle region (orb2009313)

Catalog Number	orb2009313
Category	Proteins
Description	RAB37 Peptide - middle region
Predicted Reactivity	Human
Form/Appearance	Lyophilized powder
Buffer/Preservatives	Lyophilized powder
Protein Sequence	SAKTGMNVELAFLAI AKELKYRAGHQADEPSFQIRDYVESQKKRSSCCSF
UniProt ID	Q8IWA7
Tested applications	WB
Application notes	This is a synthetic peptide designed for use in combination with RAB37 Rabbit Polyclonal Antibody (orb582864). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
NCBI	NP_783865
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)