

Product Datasheet

C4B Peptide - N-terminal region (orb2016980)

Description	C4B Peptide - N-terminal region
Tested Applications	WB
Preservatives	Lyophilized powder
Form/Appearance	Lyophilized powder
Storage	Add 100ul of sterile PBS. Final peptide concentration is 1 mg/ml in PBS. For longer periods of storage, store at -2°C. Avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	This is a synthetic peptide designed for use in combination with C4B Rabbit Polyclonal Antibody (orb578241). It may block above mentioned antibody from binding to its target protein in western blot and/or immunohistochemistry under proper experimental settings.
MW	33kDa
Uniprot ID	Q6U2E9
NCBI	NP_001002029
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)