

Product Datasheet

TRADD Rabbit Polyclonal Antibody (Biotin) (orb2081476)

Catalog Number	orb2081476
Category	Antibodies
Description	TRADD Rabbit Polyclonal Antibody (Biotin)
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Predicted Reactivity	Canine, Equine, Human
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human TRADD
Protein Sequence	Synthetic peptide located within the following region: YEQAFQLRRFVQAEGRRTLQRLVEALEENELTS LAEDLLGLTDPNGGL
UniProt ID	Q15628
MW	34kDa
Tested applications	WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

NCBI

NP_003780

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)