

Product Datasheet

CX3CR1 Rabbit Polyclonal Antibody (FITC) (orb2083527)

Catalog Number	orb2083527
Description	CX3CR1 Rabbit Polyclonal Antibody (FITC)
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human CX3CR1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EKFRRLYHLYGKCLAVLCGRSVHVDFFSSSESQSRHGSVLSSNFTYHTS
Clonality	Polyclonal
MW	39kDa
Uniprot ID	P49238

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

NCBI

NP_001328

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)