

## Product Datasheet

### NCR2 Rabbit Polyclonal Antibody (HRP) (orb2093111)

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | NCR2 Rabbit Polyclonal Antibody (HRP)                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Conjugation</b>         | HRP                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | WB                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>           | The immunogen is a synthetic peptide directed towards the C terminal region of human NCR2                                                                                                                                                                                                                                                                                                                                             |
| <b>Preservatives</b>       | Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.                                                                                                                                                                                                                                                                                                                                   |
| <b>Form/Appearance</b>     | Liquid. Purified antibody is supplied in high phosphate PBS, 100 mM phosphate, 150 mM NaCl, pH 7.6.                                                                                                                                                                                                                                                                                                                                   |
| <b>Storage</b>             | All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Sequence</b>    | Synthetic peptide located within the following region:<br>LVAKSLVLSALLVWWGDIWWKTMELRSLDTQKATCHLQQVTDLPWTSVS                                                                                                                                                                                                                                                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>MW</b>                  | 31 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Uniprot ID</b>          | <b>O95944</b>                                                                                                                                                                                                                                                                                                                                                                                                                         |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

NCBI

NP\_004819

Expiration Date

12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)