

Product Datasheet

PDIA3 Rabbit Polyclonal Antibody (FITC) (orb2093487)

Catalog Number	orb2093487
Category	Antibodies
Description	PDIA3 Rabbit Polyclonal Antibody (FITC)
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	FITC
Predicted Reactivity	Bovine, Canine, Equine, Goat, Guinea pig, Human, Mouse, Rabbit, Rat, Sheep, Zebrafish
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Purification	Affinity Purified
Protein Sequence	Synthetic peptide located within the following region: YFSPANKKLNPKKYEGGRELSDFISYLRQREATNPPVIQEEKPKKKKKAQE
UniProt ID	P30101
MW	56kDa
Tested applications	IHC, IP, WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

NCBI

NP_005304

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)