

Product Datasheet

Traf3ip1 Rabbit Polyclonal Antibody (Biotin) (orb2105236)

Catalog Number	orb2105236
Category	Antibodies
Description	Traf3ip1 Rabbit Polyclonal Antibody (Biotin) is a Biotin, affinity purified antibody which is suitable for WB application. The immunogen is a synthetic peptide directed towards the N-terminal region of Mouse Traf3ip1. This antibody is predicted to react with Bovine, Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish. It is supplied as a liquid in 1x PBS buffer.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Predicted Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Purification	Affinity Purified
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Mouse Traf3ip1
Protein Sequence	Synthetic peptide located within the following region: VIRITGFMKGLYTDAEMKSENVKDKDAKISFLQKAIDVMMVSGEPLAAK
UniProt ID	Q149C2
MW	69kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Tested applications

WB

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

NCBI**NP_082994****Expiration Date**

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.comPhone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)**Biorbyt LLC**

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.comPhone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)