

Product Datasheet

HSPB8 Rabbit Polyclonal Antibody (FITC) (orb2105652)

| | |
|-----------------------------|---|
| Catalog Number | orb2105652 |
| Category | Antibodies |
| Description | HSPB8 Rabbit Polyclonal Antibody (FITC) |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Conjugation | FITC |
| Predicted Reactivity | Bovine, Canine, Equine, Goat, Guinea pig, Human, Mouse, Rabbit, Rat, Zebrafish |
| Form/Appearance | Liquid. Purified antibody supplied in 1x PBS buffer. |
| Buffer/Preservatives | Liquid. Purified antibody supplied in 1x PBS buffer. |
| Purification | Affinity Purified |
| Immunogen | The immunogen is a synthetic peptide directed towards the middle region of human HSPB8 |
| Protein Sequence | Synthetic peptide located within the following region: PEELMVKTKDGYVEVSGKHEEKQQEGGIVSKNFTKKIQLPAEVDPTVFA |
| UniProt ID | Q9UJY1 |
| MW | 21 |
| Tested applications | WB |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

NCBI

NP_055180

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)