

## Product Datasheet

### FBXL4 Rabbit Polyclonal Antibody (FITC) (orb2105967)

|                             |                                                                                                            |
|-----------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2105967                                                                                                 |
| <b>Category</b>             | Antibodies                                                                                                 |
| <b>Description</b>          | FBXL4 Rabbit Polyclonal Antibody (FITC)                                                                    |
| <b>Clonality</b>            | Polyclonal                                                                                                 |
| <b>Species/Host</b>         | Rabbit                                                                                                     |
| <b>Conjugation</b>          | FITC                                                                                                       |
| <b>Predicted Reactivity</b> | Bovine, Canine, Equine, Human, Mouse, Porcine, Rabbit, Rat                                                 |
| <b>Form/Appearance</b>      | Liquid. Purified antibody supplied in 1x PBS buffer.                                                       |
| <b>Buffer/Preservatives</b> | Liquid. Purified antibody supplied in 1x PBS buffer.                                                       |
| <b>Purification</b>         | Affinity Purified                                                                                          |
| <b>Immunogen</b>            | The immunogen is a synthetic peptide directed towards the N terminal region of human FBXL4                 |
| <b>Protein Sequence</b>     | Synthetic peptide located within the following region:<br>VLETYHPGAVIRILACSANPYSNPPEVRWEILWSERPTKVNASQARQF |
| <b>UniProt ID</b>           | <b>Q9UKA2</b>                                                                                              |
| <b>MW</b>                   | 70kDa                                                                                                      |
| <b>Tested applications</b>  | WB                                                                                                         |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Storage**

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

**Note**

For research use only

**NCBI**

**NP\_036292**

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)