

Product Datasheet

MMP3 Rabbit Polyclonal Antibody (FITC) (orb2122800)

Catalog Number	orb2122800
Category	Antibodies
Description	MMP3 Rabbit Polyclonal Antibody (FITC) is an FITC conjugated, affinity purified antibody which is suitable for IHC, WB applications. The immunogen is a synthetic peptide directed towards the middle region of human MMP3. This antibody is predicted to react with Bovine, Canine, Equine, Human, Mouse, Porcine, Rabbit, Rat. It is supplied as a liquid in 1x PBS buffer.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	FITC
Predicted Reactivity	Bovine, Canine, Equine, Human, Mouse, Porcine, Rabbit, Rat
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Purification	Affinity Purified
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human MMP3
Protein Sequence	Synthetic peptide located within the following region: AEDFPGIDSKIDAVFEEFGFFYFTGSSQLEFDPNAKKVTHTLKSNSWLN
UniProt ID	P08254
MW	43kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Tested applications	IHC, WB
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
NCBI	NP_002413
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)