

Product Datasheet

MCM4 Rabbit Polyclonal Antibody (Biotin) (orb2131532)

Catalog Number	orb2131532
Category	Antibodies
Description	MCM4 Rabbit Polyclonal Antibody (Biotin) is a Biotin conjugated, antibody which is suitable for IHC, WB applications. The immunogen is a synthetic peptide directed towards the middle region of human MCM4. This antibody is predicted to react with Bovine, Canine, Equine, Guinea pig, Human, Mouse, Rabbit, Rat, Zebrafish. It is supplied as a liquid in 1x PBS buffer.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Biotin
Predicted Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Mouse, Rabbit, Rat, Zebrafish
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Buffer/Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human MCM4
Protein Sequence	Synthetic peptide located within the following region: VYKTHIDVIHYRKTDAAKRLHGLDEEAEQKLFSEKRVLLKELSRKPDIE
UniProt ID	P33991
MW	95kDa
Tested applications	IHC, WB

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4°C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

Note

For research use only

NCBI

NP_005905

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)