

Product Datasheet

Human ADAMTS4 ELISA Kit (orb215872)

| | |
|----------------------------|---|
| Catalog Number | orb215872 |
| Category | Assays and Kits |
| Description | Human ADAMTS4 ELISA Kit |
| Target | A disintegrin and metalloproteinase with thrombospondin motifs 4 |
| Reactivity | Human |
| Immunogen | Expression system for standard: CHO; Immunogen sequence: F213-C685 |
| UniProt ID | 075173 |
| MW | 90197 Da |
| Tested applications | ELISA |
| Specificity | Natural and recombinant human ADAMTS4 |
| Cross Reactivity | There is no detectable cross-reactivity with other relevant proteins. |
| Range | 625 pg/ml - 40,000 pg/ml |
| Sensitivity | 20 pg/ml |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Human ADAMTS4 ELISA Kit

Human ADAMTS4 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)