

## Product Datasheet

### Goat Anti-Rabbit IgG (H&L) - HRP (orb216278)

|                            |                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb216278                                                                                                       |
| <b>Description</b>         | Goat polyclonal antibody against Rabbit IgG (H&L) conjugated to HRP                                             |
| <b>Species/Host</b>        | Goat                                                                                                            |
| <b>Reactivity</b>          | Rabbit                                                                                                          |
| <b>Conjugation</b>         | HRP                                                                                                             |
| <b>Tested Applications</b> | ELISA, IH, WB                                                                                                   |
| <b>Immunogen</b>           | Rabbit IgG                                                                                                      |
| <b>Preservatives</b>       | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin          |
| <b>Form/Appearance</b>     | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                |
| <b>Note</b>                | For research use only                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                      |
| <b>Antibody Type</b>       | Secondary Antibody                                                                                              |
| <b>Source</b>              | Goat                                                                                                            |
| <b>Dilution Range</b>      | ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                 |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)