

Product Datasheet

Rabbit Anti-Goat IgG (H&L) - Biotin (orb216326)

Description	Rabbit polyclonal antibody against Goat IgG (H&L) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Biotin
Tested Applications	ELISA, IH, WB
Immunogen	Goat IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Antibody Type	Secondary Antibody
Source	Rabbit
Dilution Range	ELISA: 1:5000-1:20000, WB: 1:5000-1:20000, IHC-P: 1:500-1:1000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)