

Product Datasheet

Rabbit Anti-Sheep IgG (H&L) Antibody (FITC) (orb216377)

Catalog Number	orb216377
Category	Antibodies
Description	The Rabbit Sheep IgG (H&L) - FITC Antibody is suitable for ELISA, FC, IF, IHC. It is a Polyclonal, FITC conjugated antibody which raised against Sheep IgG. Purification: The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	FITC
Reactivity	Sheep
Form/Appearance	Liquid
Buffer/Preservatives	0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Sheep IgG
Tested applications	ELISA, FC, IF, IHC
Dilution range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Sheep IgG. No antibody was detected against non immunoglobulin serum proteins.
Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)