

Product Datasheet

ACTG1 peptide (orb216873)

Catalog Number	orb216873
Category	Proteins
Description	The peptide is used to block Anti-Gamma-actin-1 Antibody (orb213517) reactivity.
Target	ACTG1
Reactivity	Bovine, Gallus, Human, Mouse, Porcine, Rabbit, Rat, Sheep
Form/Appearance	Lyophilized powder
Purity	>85%
UniProt ID	P63260, P60709, P60710, P63259, P63261, P60711, Q9BYX7
Tested applications	Blocking
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Source	Synthetic
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Entrez	11465, 71, 100361457, 60, 81822, 11461
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)