

Product Datasheet

NRAS peptide (orb217412)

Catalog Number	orb217412
Category	Proteins
Description	The peptide is used to block Anti-NRAS/HRAS/KRAS Antibody (orb214056) reactivity.
Reactivity	Human, Mouse, Rat
Form/Appearance	Lyophilized powder
Purity	>85%
UniProt ID	P32883, P20171, P08644, P08556, P01116, Q61411, Q04970, P01111, P01112
Tested applications	Blocking
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Source	Synthetic
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Entrez	4893, 24605, 293621, 24525, 16653, 3845, 15461, 3265

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)