

Product Datasheet

Goat Cat IgG IgA IgM (heavy and light chains) Antibody (orb21757)

Catalog Number	orb21757
Category	Antibodies
Description	Precipitating polyclonal goat antiserum to cat immunoglobulins, heavy and light chains. Human standard Serum with assigned values for IgG subclasses. ... The determination of the concentration of each of the four subclasses was carried out by single radial immunodiffusion (RID, Mancini) in agar gel using subclass-specific polyclonal antisera and the WHO reference serum 67/97 as the primary standard. The nephelometric determination of the concentration of each of the four subclasses was carried out on a Beckman Coulter IMAGE
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Unconjugated
Reactivity	Feline
Form/Appearance	Delipidated, heat inactivated, lyophilized stable whole serum No preservative added, as it may interfere with the antibody activity. No foreign protein added. Total protein and IgG concentration in the antiserum are comparable to those of pooled goat seru
RRID	AB_10753859
Specificity	Precipitating polyclonal Goat antiSerum to Cat immunoglobulins, heavy and light chains
Antibody Type	Secondary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Source	Purified normal IgG and fractions containing IgA and IgM and isolated from cat serum. Freund's complete adjuvant is used in the first step of the immunization procedure.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)