

Product Datasheet

SFTPB peptide (orb217924)

Catalog Number	orb217924
Category	Proteins
Description	The peptide is used to block Anti-Surfactant Protein B Antibody (orb214568) reactivity.
Target	SFTPB
Reactivity	Human
Form/Appearance	Lyophilized powder
Purity	>85%
UniProt ID	P07988
Tested applications	Blocking
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Source	Synthetic
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Entrez	6439
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558