

Product Datasheet

ADAM17 peptide (orb217996)

| | |
|----------------------------|---|
| Catalog Number | orb217996 |
| Category | Proteins |
| Description | The peptide is used to block Anti-CD156b Antibody (orb214640) reactivity. |
| Target | ADAM17 |
| Reactivity | Human, Mouse, Rat |
| Form/Appearance | Lyophilized powder |
| Purity | >85% |
| UniProt ID | Q9Z1K9, P78536, Q9Z0F8 |
| Tested applications | Blocking |
| Application notes | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Source | Synthetic |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Entrez | 11491, 57027, 6868 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)