

Product Datasheet

UBE2A peptide (orb218066)

Catalog Number	orb218066
Category	Proteins
Description	The peptide is used to block Anti-UBE2A/B Antibody (orb214710) reactivity.
Target	UBE2A; UBE2B
Reactivity	Gallus, Human, Mouse, Rat, Zebrafish
Form/Appearance	Lyophilized powder
Purity	>85%
UniProt ID	P49459, P63147, P63146, Q9Z255, P63149
Tested applications	Blocking
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Source	Synthetic
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Entrez	22210, 81816, 22209, 7320
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)