

## Product Datasheet

### Mouse CD354 ELISA Kit (orb219637)

<b>Catalog Number</b>	orb219637
<b>Category</b>	Assays and Kits
<b>Description</b>	Sandwich ELISA for quantitative detection of Mouse CD354 concentrations.
<b>Target</b>	TREM1
<b>Reactivity</b>	Mouse
<b>UniProt ID</b>	<b>Q9JKE2</b>
<b>Tested applications</b>	ELISA
<b>Application notes</b>	Detection Method: Colorimetric
<b>Assay Type</b>	Detection
<b>Range</b>	31.2 pg/ml - 2000 pg/ml
<b>Sensitivity</b>	2 pg/ml
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Entrez</b>	<b>58217</b>
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Representative Standard Curve

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)