

Product Datasheet

Rabbit Mouse IgG1 IgG2a IgG2b IgG3 IgA IgM (heavy and light chains) (non Human), conjugated with Horseradish peroxidase Antibody (orb22093)

Catalog Number orb22093

Category Antibodies

Description Horseradish peroxidase-conjugated IgG fraction of polyclonal rabbit antiserum to

Mouse IgG, IgA and IgM, heavy and light chains, adsorbed for human cross-reactivity. In enzyme-immunocytochemical and immunohistochemical staining of monoclonal mouse monoclonal antibodies in the presence of human cell or serum components at the cellular and subcellular level. This immunoconjugate is not pre-diluted. The optimum working dilution of each conjugate should be established by titration before being used. Excess labelled antibody must be avoided because it may cause high unspecific background staining and interfere with the specific signal. Working dilutions for histochemical and cytochemical use are usually between 1:100 and 1:250; in ELISA and comparable non-

precipitating antibody-binding assays between 1:500 and 1:5,000.

Clonality Polyclonal

Species/Host Rabbit

Conjugation HRP

Reactivity Mouse

Form/Appearance Peroxidase-coupled purified hyperimmune rabbit IgG lyophilized from a solution

in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may

interfere with the antibody activity.

RRID AB_10932451

Antibody Type Secondary Antibody

Phone: <u>+1 (415) 906-5211</u> | Fax: <u>+1 (415) 651-8558</u>





Source

Purified polyclonal mouse IgG, and homogenous IgA and IgM isolated from mouse serum. Freund's complete adjuvant is used in the first step of the immunization procedure.

Storage

The lyophilized conjugate is shipped at ambient temperature and may be stored at $+4^{\circ}$ C; prolonged storage at or below -20° C. Spun down to remove insoluble particles, divided into small aliquots, frozen and stored at or below -20° C. Prior to use, an aliquot is thawed slowly in the dark at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at $+4^{\circ}$ C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the immunoconjugate.

Note For research use only

Expiration Date 12 months from date of receipt.