

Product Datasheet

CRISPR/Cas9 Monoclonal Antibody (orb223765)

Catalog Number	orb223765
Category	Antibodies
Description	CRISPR/Cas9 Monoclonal Antibody
Target	CRISPR/Cas9
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1k
Conjugation	Unconjugated
Reactivity	Bacteria
Concentration	1mg/ml
Buffer/Preservatives	Purified antibody (from supernatant) containing PBS + 0.09% sodium azide
Purification	Protein A
Immunogen	This antibody was raised against a recombinant protein within the N-terminal region of Streptococcus pyogenes Cas9
UniProt ID	q99zw2
Tested applications	ICC, WB
Dilution range	Optimal antibody dilution should be determined by titration, however as a guideline try at; IB 0.1-1µg/ml
Application notes	X63.Ag8-653

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Clone Number	7A9-3A3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Entrez	69900935
Expiration Date	12 months from date of receipt.

anti S. Pyogenes CRISPR/Cas9, clone 7A9

Western blot analysis of using Caspase 9 antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)