

Product Datasheet

HDAC2 Monoclonal Antibody [3F3] - Fluorescein (FITC) / 50ul (orb223781)

Description	Mouse monoclonal antibody against HDAC2 conjugated to FITC
Species/Host	Mouse
Reactivity	Bovine, Human, Rat
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF
Immunogen	Synthetic peptide, corresponding to amino acids 473-488 of Human HDAC2
Target	HDAC2
Preservatives	Purified antibody conjugated to Fluorescein using Innova Biosciences Lightning-Link™, supplied in Phosphate buffered saline (PBS) containing 0.09% sodium azide
Concentration	1 mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	3F3
Uniprot ID	Q92769, P56519, P70288
Entrez	395635, 84577, 15182, 407223, 3066
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)