

Product Datasheet

DALDA TFA (orb2277395)

| | |
|------------------------|--|
| Catalog Number | orb2277395 |
| Category | Proteins |
| Description | DALDA TFA is a highly selective and potent μ -opioid receptor agonist with a K_i of 1.69 nM. It is used in research to study opioid receptor pharmacology, including investigations of antinociception and respiratory depression in both in vitro binding assays and in vivo animal models. |
| Target | Opioid Receptor |
| MW | 611.74 (free base) |
| Formula | $C_{30}H_{45}N_9O_5 \cdot xC_2HF_3O_2$ |
| Storage | -20°C |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)