

## Product Datasheet

### RORyt inverse agonist 31 (orb2277919)

<b>Catalog Number</b>	orb2277919
<b>Category</b>	Small Molecules
<b>Description</b>	RORyt inverse agonist 31 (14g) is a potent antagonist of the retinoic acid receptor-related orphan receptor $\gamma$ t (RORyt), with an inhibitory concentration (IC <sub>50</sub> ) of 0.428 $\mu$ M. It has shown efficacy in reducing Imiquimod-induced psoriasis severity in murine models.
<b>Target</b>	ROR
<b>MW</b>	555.36
<b>Formula</b>	$C_{23}H_{15}Cl_2F_3N_4O_3S$
<b>Storage</b>	-20°C
<b>Note</b>	For research use only

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)