

Product Datasheet

A3AR agonist 2 (orb2277994)

Catalog Number	orb2277994
Category	Small Molecules
Description	Compound 19, a selective A3 Adenosine Receptor (A3AR) agonist (K _i : 22.1 nM), effectively stimulates β-arrestin2 recruitment with an EC ₅₀ value of 4.36 nM. This compound is utilized in research pertaining to inflammatory diseases, ischemia, cancer, neuropathic pain, liver diseases, and other chronic conditions.
Target	Adenosine Receptor
MW	622.67
Formula	C ₃₄ H ₃₄ N ₆ O ₆
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)