

Product Datasheet

Blinatumomab Antibody (orb2278710)

Catalog Number	orb2278710
Category	Antibodies
Description	Blinatumomab is a bispecific T-cell engager antibody that targets CD19 on B-cells and CD3 on T-cells, redirecting T-cell cytotoxicity to malignant cells. It is utilized in research models of B-cell malignancies, including acute lymphoblastic leukemia and non-Hodgkin's lymphoma, for both in vitro cytotoxicity assays and in vivo efficacy studies.
Target	Integrin,Immunology/Inflammation related
Clonality	Recombinant
Conjugation	Unconjugated
Concentration	4.29 mg/mL
Purity	>90% (May vary between batches)
MW	54 kDa
Antibody Type	Inhibitor Antibody
Target Areas	Integrin Immunology/Inflammation related
CAS Number	853426-35-4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558