

Product Datasheet

Isocitrate Dehydrogenase 1 (IDH1) Mouse Monoclonal Antibody (Capture) (orb2280071)

Catalog Number	orb2280071
Category	Antibodies
Description	Isocitrate Dehydrogenase 1 (IDH1) Mouse Monoclonal Antibody (Capture)
Clonality	Monoclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Purified,Liquid
Buffer/Preservatives	PBS, pH 7.4
Purity	>95% pure (SDS-PAGE)
Purification	Protein A Chromatography
MW	150 kDa
Tested applications	CLIA
Specificity	Specific for IDH1
Matched Pairs	Capture: orb2280071; Detection: orb2280072;
Source	Mammalian cell
Sample Types	Test sample type: Serum / Plasma. Test sample source: Human

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)