

Product Datasheet

Isocitrate Dehydrogenase 1 (IDH1) Mouse Monoclonal Antibody (Capture) (orb2280071)

| | |
|-----------------------------|---|
| Catalog Number | orb2280071 |
| Category | Antibodies |
| Description | Isocitrate Dehydrogenase 1 (IDH1) Mouse Monoclonal Antibody (Capture) |
| Clonality | Monoclonal |
| Isotype | IgG |
| Conjugation | Unconjugated |
| Reactivity | Mouse |
| Form/Appearance | Purified,Liquid |
| Buffer/Preservatives | PBS, pH 7.4 |
| Purity | >95% pure (SDS-PAGE) |
| Purification | Protein A Chromatography |
| MW | 150 kDa |
| Tested applications | CLIA |
| Specificity | Specific for IDH1 |
| Matched Pairs | Capture: orb2280071; Detection: orb2280072; |
| Source | Mammalian cell |
| Sample Types | Test sample type: Serum / Plasma. Test sample source: Human |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)