


Product Datasheet

Human ARHGAP1 (Rho GTPase-activating protein 1) ELISA Kit (orb2280369)

Description	Human ARHGAP1 (Rho GTPase-activating protein 1) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	ARHGAP1
Storage	2 - 8 °C for 12 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Assay Time	4 hours
Uniprot ID	Q07960
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)