

Product Datasheet

BR Prestained Protein Ladder B (10-180 kDa) (orb2280629)

| | |
|------------------------|---|
| Catalog Number | orb2280629 |
| Category | Proteins |
| Description | <p>The BR Prestained Protein Ladder B, is a combination of 10 pre-stained proteins with molecular weights from 11 to 180 kDa. The 10 recombinant proteins are covalently coupled with blue chromophore, while 1 green band at 25 kDa and a red band at 75 kDa, serve as reference bands. The Prestained Protein Ladder B, keeps track of the size and separation of proteins during SDS-polyacrylamide gel electrophoresis, approximating protein size, and validating Western transfer efficiency on PVDF, nylon, or nitrocellulose membranes.</p> |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)