

Product Datasheet

IRIS9 Prestained Protein Ladder (9 pre-stained bands, 15-180 kDa) (orb2280632)

Catalog Number orb2280632

Description The BR Prestained Protein Ladder E is a combination of 9 pre-stained proteins with molecular weights from 15 to 180 kDa. The 9 recombinant proteins are covalently coupled with blue chromophore, while 2 red bands at 70 kDa and 30 kDa, a green band at 15 kDa and a newly designed peacock green band at 60 kDa serve as reference bands. The BR Prestained Protein Ladder E keeps track of the size and separation of proteins during SDS-polyacrylamide gel electrophoresis, approximating the target protein size and validating the Western transfer efficiency on PVDF, nylon, or nitrocellulose membranes.

Note For research use only

Expiration Date 6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)