

Product Datasheet

MOUSE RABBIT CD45:RPE Antibody (orb2281278)

Catalog Number	orb2281278
Category	Antibodies
Description	Mouse anti Rabbit CD45:RPE
Species/Host	Mouse
Isotype	IgG1
Conjugation	RPE
Reactivity	Rabbit
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose; Buffer: Phosphate buffered saline
Purification	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Immunogen	Glycoproteins isolated from the T cell line, RL-5.
Tested applications	FC
Dilution range	Min Dilution: Neat, Max Dilution: 1/10. Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul
Specificity	CD45
Antibody Type	Primary Antibody
Clone Number	L12/201

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)