

Product Datasheet

Rabbit Anti-Human IgG H&L, BF405 conjugated (orb2283904)

Catalog Number	orb2283904
Category	Antibodies
Description	The Rabbit Anti-Human IgG H&L, BF405 conjugated is a BF405 conjugated polyclonal antibody suitable for FC, ICC, IF. This product reacts with Human samples.
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	BF405
Reactivity	Human
Purification	Affinity purified by Protein A
Immunogen	Native Human IgG full length
Tested applications	FC, ICC, IF
Dilution range	ICC/IF=1:100-1000, IF=1:100-1000, Flow-Cyt=1:100-1000
Antibody Type	Secondary Antibody
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)